


Choose how to wear your help button.

Decide if you want to wear your Help Button as a wristband or pendant.

Next, snap the button into the wristband or clip the necklace onto the top of the button.



Help Button with Fall Detection *(optional)*



FALL DETECTION IS AN OPTIONAL SERVICE THAT CAN BE PROVIDED WITH YOUR SYSTEM. PLEASE CALL 800.998.2400 FOR MORE INFORMATION.


You should press the HELP BUTTON if you need help. If you are unable to press the button and a fall is detected by the HELP BUTTON, the system beeps for 15 seconds and then sends an emergency fall alert to the monitoring center.

- To CANCEL a false fall detection alert, immediately press and hold the HELP BUTTON for 8 seconds.
- For best results, wear your HELP BUTTON **OUTSIDE** of all your clothes.
- Handle your HELP BUTTON with care when putting it on or taking it off, as you may accidentally activate it.


Fall Detection does not detect 100% of falls. If able, users should always push their button when they need help.

Getting to know your console.

If the **GREEN POWER** light is **NOT ON** please call Customer Service.



If the **YELLOW TROUBLE** light **IS ON** please call Customer Service.



Your Help Button is waterproof and **CAN** be worn in the shower.



Press the **CLEAR** button to temporarily silence a voice announcement.



Protect your medical alert equipment.




Lost or misplaced help button or base station



Damage due to dust, heat and humidity



Drops, spills or accidents



Power surges

Our medical alert Assurance Plan provides you with additional peace of mind knowing your lifesaving medical alert equipment is always protected from everyday mishaps. Lose your help button? Grandchildren spill soda on the medical alert equipment? Don't worry, you are protected with our Assurance Plan.

Call 800.998.2400 for more information or to enroll, now.

IMPORTANT:

Sign the Service Agreement where noted and return it in the enclosed postage-paid envelope.

LF 6535 03/20 ©LifeStation 2020



If you have any questions while setting up your system, call (800) 998-2400 for assistance.

If you have any questions while setting up your system, call (800) 998-2400 for assistance.

If you have any questions while setting up your system, call (800) 998-2400 for assistance.

Getting Started

- 1 Follow the steps inside to setup your system.
- 2 Review your enclosed **Account Information Form** and call us if you need to make any changes.
- 3 Sign the **Service Agreement** where noted and return it in the enclosed postage-paid envelope.


Start here.

First, determine where your console will be placed.

Your console needs to be near an **electrical outlet that is not controlled by a switch** and near a **telephone** connected to a wall jack.

Next, bring the box with **all of the contents** to this location. You are now ready to set up your console.

Now, use the following steps to the right in order to connect and activate your medical alert system.



If you have any questions while setting up your system, call **(800) 998-2400** for assistance.

IMPORTANT: Sign the Service Agreement where noted and return it in the enclosed postage-paid envelope.